

## Just Point and See:

Near and distance viewing is made easy at school, work, or home. Everything is within your reach with Flipper's innovative rotating camera. Connect Flipper to any TV or monitor, then roll and point at any object. With Flipper enjoy full-color auto focus magnification and multiple contrast options. With so many easy to use features Flipper is the ultimate in flexibility.


### Flipper Features:

- Compact and fully portable
- Lightweight at just over 1 lb.
- True Color with auto focus
- 225 degree camera rotation to view images at any distance
- Distance, intermediate and near viewing
- 6 to 50x adjustable magnification
- Combine Flipper with a 12" battery operated LCD or connect to any TV
- Computer and laptop compatible (additional hardware required)
- 4 viewing modes to optimize contrast and brightness
- Combine Flipper with an optional stand to create a desktop magnifier



“ Today I am able to work again with Flipper. The rotating camera enables me to perform daily tasks more effectively. By rotating the camera head, I can read, write and view the titles on the library aisles.”

–Damien Utkewicz  
Cataloger, Princeton Theological Seminary  
Princeton, NJ



*Flipper Panel*

**enhanced**  
**vision**

**(888) 811-3161**

**(714) 465-3400**

[www.enhancedvision.com](http://www.enhancedvision.com)



## Flipper & USB Bridge

Pair Flipper with the USB Bridge and experience portability like never before. Use your laptop as a viewing screen.



*Flipper Stand*

## Flipper Stand

Use Flipper with the optional Stand and use your own TV to create a desktop system.

To schedule a no obligation demonstration, please call (888) 811-3161.

Local Dealer:



*Flipper Panel*

## Flipper Panel

The Flipper Panel incorporates a 12" LCD screen, with the Flipper attached to a rechargeable battery pack. The battery pack offers over 5 hours of extended use. View full-color, magnified images on a slim, portable screen.

*Flipper Panel  
Backpack Carrying  
Case included*



**enhanced**  
**vision**

**(888) 811-3161**

**(714) 465-3400**

[www.enhancedvision.com](http://www.enhancedvision.com)